



*Halogen lamp*

**174 range**

**174-BLANC**

**CARINA**  
*Medical*



**CE**

ISO 9001  
ISO 13485

## TECHNICAL DETAILS

- Halogen lamp 50 Watts, white colour
- Flexible arm, length 650mm
- Light intensity: 870 Lumens - 3000° Kelvin (neutral white)
- CRI: 100/100
- Luminous power: 40000 Lux/ 50cm, illuminated diameter: 10cm
- Power cable length 3 meters
- Class I medical device

## OPTIONS

- No contact switch
- Mobile base with wheels ( ACCI73-PIED)

Powerful and versatile,  
this lamp is ideal for  
general medicine and  
minor surgery



### WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years



CONTACT

CARINA  
*Médical*

Parc d'entreprises Visionis  
278 rue de l'industrie  
01090 GUÉREINS- FRANCE

**COMMERCIAL SERVICE**

Tel : +33(0)4 74 66 69 49  
Fax : +33(0)4 74 66 69 50  
[export@carinamedical.com](mailto:export@carinamedical.com)