



Infusion stand
315 range
315-

CARINA
Médical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Folding stainless steel infusion stand, 2 polyamid hooks
- Base diameter: 55cm
- Adjustable height from 1.05 to 1.5m
- Dimension when folded: 1.10m
- Supplied with transport bag

OPTIONS

- No contact switch
- Mobile base with wheels (ACCI73-PIED)

Transport bag
include



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 1 year



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com