



*Polyurethane chair
64 range*

6401 CHR

CARINA
Medical



CE

ISO 9001
ISO 13485

TECHNICAL DETAILS

- Height adjustable chair with chromed gas column
- Manual control handle under the seat
- With very comfortable ergonomic polyurethane backrest and seat
- Equipped with 5 wheels diam 50mm
- Chromed base with 5 branches, diameter 55cm

DIMENSIONS

- Adjustable height from 47 to 67cm
- Load capacity: 130kg in dynamic

You may also like :

Ref. **6411CHR**: Same as 6401CHR but with fixed skids instead of wheels



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for the upholstery)

ACCESSORIES



ACC 14
Foot control




ACC 15
Footrest




ACC 16
Chromed base, diameter 38cm




ACC ROULETTE-02
5 single -wheels, diameter 50mm




ACC ROULETTE-03
5 double-wheels for hard floors,
diameter 50mm



ACC ROULETTE-04
5 anti static wheels, diameter 50mm



ACC ROULETTE-05
5 auto-braking wheels, in load or
out of load



HH/HB
High or low height



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 2 years



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE
COMMERCIAL SERVICE

Tel : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com