




Halogen lamp

173 range
173-BLANC

CARINA
Médical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Halogen lamp 20 Watts, white colour
- Flexible arm, length 650mm
- Light intensity: 300 Lumens - 3000° Kelvin (neutral white)
- CRI: 100/100
- Luminous power: 50000 Lux/ 30cm, illuminated diameter: 5cm
- Power cable length 3 meters
- Class I medical device

OPTIONS

- No contact switch
- Mobile base with wheels (ACCI73-PIED)

Powerful and versatile,
this lamp is ideal for
general medicine and
minor surgery



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days
upon receipt of order
- ✓ Warranty: 5 years



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com