



*Picture of the model with optionnal headrest cushion ACCCOUSSIN-01
option 4 wheels with central brake system ACC150-MULT and
ACCMUL-RPIEDS Electric foot rest with cushion*

Electric multi-treatment chair

MULTIKA version V2 range

MULT-V2-75-ME2

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 3 sections multi-treatment chair, width 75cm
- Control through multifunction hand remote with buttons for memory positions (Seat position and transfer position)
- Adjustable backrest and legrest with electric jacks
- Tredelenburg position
- One section mattress , thickness 90mm (6cm of HR35 foam + 3cm of viscoelastic memory foam)
- Equipped with a pair of armrest on knee system, synchronized with the backrest, a push bar into the backrest, and a pair of DIN rails 25x10 to fix accessories
- Shapped base for ideal cleaning, base colour light grey (RAL 7035)
- On 4 design wheels with one brake on each wheel diam 125mm

DIMENSIONS

- Width 75cm, width overall 85cm
- Adjustable height from 56 to 96cm
- Load capacity: 200kg in dynamic



You may also like :

Réf. **MULT-V2-65-ME2**: Same as MULT-V2-75-ME2 but width 65cm



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years (1 year for the upholstery)

ZOOM MULT-V2-75-ME2 with accessories




ACCESSORIES



ACC MUL-BPPIEDS

Retractable push handles on the foot side



ACC 125

Paper roll holder for exam sheet



ACC BATTERIE01

Battery, capacity 20 cycles
(1 cycle = 1x height's up/ down +
1x backrest's up/down)



ACC COUSSIN-01

Comfortable cushion
50 x 40cm



ACC CLAMEDIAM16

Adaptation clamp in aluminium
for fixation on rail dim 25x10, the unit
siam.14mm, diam.16mm



ACC MULT-320INO1

Adjustable infusion stand with 2 plastic
hooks, without clamp



ACC 514-MULT75ME

Plastic protection for leg rest



ACC MULTIK08-AM

Oxygen bottle holder with fixation
straps



ACC 150-MULT

4 wheels with central brake system
diam 150mm



ACC-MUL-RPIEDS

Electric foot rest with cushion



ACC MULTIK02-PIED

Emergency pedal to be put in
Trendelenburg position quickly



ACC EQUIPOTENT01

Equipotentiality plug



ACC MUL-5ROUES

5th wheel to optimize transfers
and mobility of the chair, to be
order with ACC150-MULT option



ACC ROUE-MOTOR01

Motorized assistance



ACC MULTIK22


Manual release of the backrest

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range

Upholstery : Classic range



Arno



Garonne




Potomac




Amour




Mississippi




St Laurent



Rhône



Missouri



Seine



Amazone

Upholstery : Urban range



Abysses



Cendre



Origan




Géranium




Réglisse




Prune




Perle



Ficelle



Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole




Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 5 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com